

Product datasheet for **GA202261**

Mouse *Kcnh2* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnh2</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily H (eag-related), member 2
Symbol:	<i>Kcnh2</i>
UniProt ID:	Q35219
Synonyms:	AI326795; ERG1; LQT; Lqt2; M-erg; Merg1; merg1a; merg1b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)