

## Product datasheet for **GA202262**

### Mouse *Kcnh3* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnh3</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily H (eag-related), member 3
Symbol:	<i>Kcnh3</i>
UniProt ID:	<a href="#">Q9WVJ0</a>
Synonyms:	AU019351; C030044P22Rik; Elk2; Kv12.2; Melk2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)