


Product datasheet for **GA202277**

Mouse *Kcnk3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnk3</i> CRISPRa kit - CRISPR gene activation of mouse potassium channel, subfamily K, member 3
Symbol:	<i>Kcnk3</i>
UniProt ID:	Q35111
Synonyms:	cTBAK-1; TASK; Task-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)