

Product datasheet for **GA202278**

Mouse *Kcnk4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnk4</i> CRISPRa kit - CRISPR gene activation of mouse potassium channel, subfamily K, member 4
Symbol:	<i>Kcnk4</i>
UniProt ID:	O88454
Synonyms:	MLZ-622; Tex40; TRAAK; TRAAKt
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)