

Product datasheet for **GA202283**

Mouse *Kcnn4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnn4</i> CRISPRa kit - CRISPR gene activation of mouse potassium intermediate/small conductance calcium-activated channel, subfamily N, member 4
Symbol:	<i>Kcnn4</i>
UniProt ID:	O89109
Synonyms:	IK1; IKCA1; KCa3.1; KCA4; mIKCa1; SK4; SKCas
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)