

Product datasheet for **GA202284**

Mouse *Kcnq1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnq1</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily Q, member 1
Symbol:	<i>Kcnq1</i>
UniProt ID:	P97414
Synonyms:	AW559127; <i>Kcna9</i> ; <i>KVLQT1</i>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)