

Product datasheet for **GA202287**

Mouse *Kcns2* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcns2</i> CRISPRa kit - CRISPR gene activation of mouse K ⁺ voltage-gated channel, subfamily S, 2
Symbol:	<i>Kcns2</i>
UniProt ID:	Q35174
Synonyms:	E130006J24Rik; Kv9.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)