

Product datasheet for **GA202318**

Mouse Kifc3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kifc3 CRISPRa kit - CRISPR gene activation of mouse kinesin family member C3
Symbol:	Kifc3
Components:	GA202318G1, Kifc3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA202318G2, Kifc3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA202318G3, Kifc3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O35231</u>
Synonyms:	AI325457; BB123200
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)