

Product datasheet for GA202320

Mouse Uhmk1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Uhmk1 CRISPRa kit - CRISPR gene activation of mouse U2AF homology motif (UHM) kinase 1
Symbol:	Uhmk1
UniProt ID:	P97343
Synonyms:	4732477C12Rik; 4930500M09Rik; AA673513; AI449218; AU021979; C820018A03Rik; KIS; Kist
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)