

## Product datasheet for **GA202333**

### Mouse Klk1b16 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Klk1b16 CRISPRa kit - CRISPR gene activation of mouse kallikrein 1-related peptidase b16                                                                                                                                                                                                                         |
| Symbol:       | Klk1b16                                                                                                                                                                                                                                                                                                          |
| Components:   | GA202333G1, Klk1b16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA202333G2, Klk1b16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA202333G3, Klk1b16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>P04071</u>                                                                                                                                                                                                                                                                                                    |
| Synonyms:     | Klk16; mGk-16                                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)