

Product datasheet for **GA202363**

Mouse **Hivep3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Hivep3 CRISPRa kit - CRISPR gene activation of mouse human immunodeficiency virus type 1 enhancer binding protein 3
Symbol:	Hivep3
UniProt ID:	A2A884
Synonyms:	2900056N03Rik; A130075N07; AI848000; E030045D18Rik; KBP-1; KBP1; Krc; Rc; Schnurri-3; Shn3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)