

Product datasheet for **GA202393**

Mouse Ktn1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ktn1 CRISPRa kit - CRISPR gene activation of mouse kinectin 1
Symbol:	Ktn1
Components:	GA202393G1, Ktn1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA202393G2, Ktn1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA202393G3, Ktn1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q61595
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)