

Product datasheet for **GA202399**

Mouse Lag3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lag3 CRISPRa kit - CRISPR gene activation of mouse lymphocyte-activation gene 3
Symbol:	Lag3
UniProt ID:	Q61790
Synonyms:	CD223; LAG-3; Ly66
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)