

Product datasheet for **GA202411**

Mouse Lamp1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lamp1 CRISPRa kit - CRISPR gene activation of mouse lysosomal-associated membrane protein 1
Symbol:	Lamp1
UniProt ID:	P11438
Synonyms:	CD107a; Lamp-1; LGP-120; LGP-A; P2B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)