

Product datasheet for **GA202412**

Mouse Lamp2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lamp2 CRISPRa kit - CRISPR gene activation of mouse lysosomal-associated membrane protein 2
Symbol:	Lamp2
UniProt ID:	P17047
Synonyms:	CD107b; Lamp-2; Lamp-2a; Lamp-2b; Lamp-2c; Lamp II; LGP-B; Mac3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)