

Product datasheet for **GA202452**

Mouse Lhx3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lhx3 CRISPRa kit - CRISPR gene activation of mouse LIM homeobox protein 3
Symbol:	Lhx3
Components:	GA202452G1, Lhx3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA202452G2, Lhx3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA202452G3, Lhx3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P50481</u>
Synonyms:	Lim3; mLim-3; mLIM3; P-LIM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)