

Product datasheet for **GA202476**

Mouse Psmb8 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Psmb8 CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) subunit, beta type 8 (large multifunctional peptidase 7) |
| Symbol: | Psmb8 |
| UniProt ID: | P28063 |
| Synonyms: | Lmp-7; Lmp7 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)