

Product datasheet for GA202499

Mouse Lrrn3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lrrn3 CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat protein 3, neuronal
Symbol:	Lrrn3
UniProt ID:	<u>Q8CBC6</u>
Synonyms:	NLRR-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)