

## Product datasheet for **GA202507**

### Mouse *Ltb4r1* activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | <i>Ltb4r1</i> CRISPRa kit - CRISPR gene activation of mouse leukotriene B4 receptor 1 |
| Symbol:       | <i>Ltb4r1</i>   |
| UniProt ID:   | <a href="#">O88855</a>  |
| Synonyms:     | BLT1; BLTR; <i>Ltb4r</i> ; mBLTR  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)               |

#### Product images:



[View online >](#)