

Product datasheet for **GA202513**

Mouse Ltk activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Ltk CRISPRa kit - CRISPR gene activation of mouse leukocyte tyrosine kinase |
| Symbol: | Ltk |
| UniProt ID: | P08923 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)