

Product datasheet for **GA202525**

Mouse Tmed1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmed1 CRISPRa kit - CRISPR gene activation of mouse transmembrane p24 trafficking protein 1
Symbol:	Tmed1
UniProt ID:	Q3V009
Synonyms:	Il1rl1l; Ly84l; St2l
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)