

Product datasheet for **GA202535**

Mouse Tm4sf1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tm4sf1 CRISPRa kit - CRISPR gene activation of mouse transmembrane 4 superfamily member 1
Symbol:	Tm4sf1
UniProt ID:	Q64302
Synonyms:	L6; M3s1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)