

Product datasheet for GA202542

Mouse Mxd3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mxd3 CRISPRa kit - CRISPR gene activation of mouse Max dimerization protein 3
Symbol:	Mxd3
UniProt ID:	<u>Q80US8</u>
Synonyms:	4631412E13Rik; bHLHc13; Mad3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)