

Product datasheet for **GA202730**

Mouse **Msx3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Msx3 CRISPRa kit - CRISPR gene activation of mouse msh homeobox 3
Symbol:	Msx3
Components:	GA202730G1, Msx3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA202730G2, Msx3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA202730G3, Msx3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P70354
Synonyms:	AI323377
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)