

Product datasheet for **GA202735**

Mouse Mt3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mt3 CRISPRa kit - CRISPR gene activation of mouse metallothionein 3
Symbol:	Mt3
Components:	GA202735G1, Mt3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA202735G2, Mt3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA202735G3, Mt3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P28184</u>
Synonyms:	Mt-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)