


Product datasheet for **GA202767**

Mouse Mup3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mup3 CRISPRa kit - CRISPR gene activation of mouse major urinary protein 3
Symbol:	Mup3
Components:	GA202767G1, Mup3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA202767G2, Mup3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA202767G3, Mup3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P04939</u>
Synonyms:	Mup-3; MUP15; Mup25; MUPIII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)