

Product datasheet for **GA203002**

Mouse Oc90 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Oc90 CRISPRa kit - CRISPR gene activation of mouse otoconin 90
Symbol:	Oc90
Components:	GA203002G1, Oc90 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203002G2, Oc90 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203002G3, Oc90 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Z0L3</u>
Synonyms:	oc95; Ocn-95; PLA2L; Pla2ll
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)