

## Product datasheet for **GA203076**

### Mouse P2rx1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P2rx1 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2X, ligand-gated ion channel, 1
Symbol:	P2rx1
UniProt ID:	<a href="#">P51576</a>
Synonyms:	AI323649; BB122383; P2x; Pdc3; RP-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)