

Product datasheet for **GA203077**

Mouse P2rx4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2rx4 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2X, ligand-gated ion channel 4
Symbol:	P2rx4
Synonyms:	AI504491; AW555605; D5Erttd444e; P2X4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)