

Product datasheet for **GA203078**

Mouse P2rx7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2rx7 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2X, ligand-gated ion channel, 7
Symbol:	P2rx7
UniProt ID:	Q9Z1M0
Synonyms:	A1467586; P2X(7); P2X7R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

