

Product datasheet for **GA203079**

Mouse P2rx6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2rx6 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2X, ligand-gated ion channel, 6
Symbol:	P2rx6
UniProt ID:	<u>O54803</u>
Synonyms:	P2rxl1; P2x6; P2xm
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)