

Product datasheet for **GA203081**

Mouse P2ry2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2ry2 CRISPRa kit - CRISPR gene activation of mouse purinergic receptor P2Y, G-protein coupled 2
Symbol:	P2ry2
UniProt ID:	P35383
Synonyms:	P2U1; P2Y2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)