

Product datasheet for **GA203090**

Mouse Pafah1b3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pafah1b3 CRISPRa kit - CRISPR gene activation of mouse platelet-activating factor acetylhydrolase, isoform 1b, subunit 3
Symbol:	Pafah1b3
UniProt ID:	Q61205
Synonyms:	mus[g]; Pafahg
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)