

## Product datasheet for **GA203097**

### Mouse Reg3b activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Reg3b CRISPRa kit - CRISPR gene activation of mouse regenerating islet-derived 3 beta
Symbol:	Reg3b
UniProt ID:	<a href="#">P35230</a>
Synonyms:	HIP; Pap; PAP1; REG-III
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)