

Product datasheet for **GA203101**

Mouse Pax3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pax3 CRISPRa kit - CRISPR gene activation of mouse paired box 3
Symbol:	Pax3
UniProt ID:	P24610
Synonyms:	Pax-3; Sp; Splchl2; splotch
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)