

## Product datasheet for **GA203102**

### Mouse Pax4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Pax4 CRISPRa kit - CRISPR gene activation of mouse paired box 4
Symbol:	Pax4
UniProt ID:	<a href="#">P32115</a>
Synonyms:	Pax-4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)