


Product datasheet for **GA203109**

Mouse Pbx2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pbx2 CRISPRa kit - CRISPR gene activation of mouse pre B cell leukemia homeobox 2
Symbol:	Pbx2
Components:	GA203109G1, Pbx2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203109G2, Pbx2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203109G3, Pbx2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	O35984
Synonyms:	AU043397; G17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)