


Product datasheet for **GA203117**

Mouse Pcm1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcm1 CRISPRa kit - CRISPR gene activation of mouse pericentriolar material 1
Symbol:	Pcm1
Components:	GA203117G1, Pcm1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203117G2, Pcm1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203117G3, Pcm1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9R0L6
Synonyms:	2600002H09Rik; 9430077F19Rik; C030044G17Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)