

Product datasheet for **GA203120**

Mouse Pcnt activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcnt CRISPRa kit – CRISPR gene activation of mouse pericentrin (kendrin)
Symbol:	Pcnt
Components:	GA203120G1, Pcnt gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203120G2, Pcnt gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203120G3, Pcnt gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P48725</u>
Synonyms:	AW476095; C86676; KEN; kendrin; m239Asp; m275Asp; Pcnt2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

