

Product datasheet for **GA203127**

Mouse Pcsk4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcsk4 CRISPRa kit - CRISPR gene activation of mouse proprotein convertase subtilisin/kexin type 4
Symbol:	Pcsk4
Components:	GA203127G1, Pcsk4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203127G2, Pcsk4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203127G3, Pcsk4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P29121</u>
Synonyms:	AI647044; L21221; NEC 3; NEC3; PC4; SPC5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)