

Product datasheet for **GA203130**

Mouse *Pcsk7* activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | <i>Pcsk7</i> CRISPRa kit - CRISPR gene activation of mouse proprotein convertase subtilisin/kexin type 7 |
| Symbol: | <i>Pcsk7</i> |
| UniProt ID: | Q61139 |
| Synonyms: | AA959856; PC7; SPC7 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)