

## Product datasheet for **GA203140**

### Mouse **Pdcd11** activation kit by **CRISPRa**

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | Pdcd11 CRISPRa kit - CRISPR gene activation of mouse programmed cell death 11                                                                                                                                                                                                                                 |
| <b>Symbol:</b>       | Pdcd11                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA203140G1, Pdcd11 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA203140G2, Pdcd11 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA203140G3, Pdcd11 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q6NS46</a></u>                                                                                                                                                                                                                                                                                 |
| <b>Synonyms:</b>     | 1110021122Rik; ALG-4; mKIAA0185; Pdcd7                                                                                                                                                                                                                                                                        |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

