

Product datasheet for **GA203144**

Mouse Pde3b activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Pde3b CRISPRa kit - CRISPR gene activation of mouse phosphodiesterase 3B, cGMP-inhibited |
| Symbol: | Pde3b |
| UniProt ID: | Q61409 |
| Synonyms: | 9830102A01Rik; AI847709 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:

