

## Product datasheet for **GA203160**

### Mouse Padi3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Padi3 CRISPRa kit - CRISPR gene activation of mouse peptidyl arginine deiminase, type III
Symbol:	Padi3
UniProt ID:	<a href="#">Q9Z184</a>
Synonyms:	AI507004; mKIAA4171; Pad3; Pdi3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)