

Product datasheet for **GA203166**

Mouse Pdx1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Pdx1 CRISPRa kit - CRISPR gene activation of mouse pancreatic and duodenal homeobox 1 |
| Symbol: | Pdx1 |
| UniProt ID: | P52946 |
| Synonyms: | IDX-1; IPF-1; Ipf1; Mody4; pdx-1; STF-1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)