

Product datasheet for GA203168

Mouse Pea15a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pea15a CRISPRa kit – CRISPR gene activation of mouse phosphoprotein enriched in astrocytes I5A
Symbol:	Pea15a
Components:	GA203168G1, Pea15a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA203168G2, Pea15a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA203168G3, Pea15a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q62048</u>
Synonyms:	Matl; PEA-15; Pea15; Pkcs15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

