

Product datasheet for **GA203171**

Mouse Peg3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Peg3 CRISPRa kit - CRISPR gene activation of mouse paternally expressed 3
Symbol:	Peg3
UniProt ID:	<u>Q3URU2</u>
Synonyms:	AL022617; ASF-1; End4; Gcap4; mKIAA0287; Pw1; Zfp102
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)