

Product datasheet for **GA203188**

Mouse Pfk1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pfk1 CRISPRa kit - CRISPR gene activation of mouse phosphofruktokinase, liver, B-type
Symbol:	Pfk1
UniProt ID:	<u>P12382</u>
Synonyms:	AA407869; ATP-PFK; PFK-L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)