

Product datasheet for **GA203197**

Mouse P_{gk2} activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P _{gk2} CRISPRa kit - CRISPR gene activation of mouse phosphoglycerate kinase 2
Symbol:	P _{gk2}
UniProt ID:	P09041
Synonyms:	P _{gk} -2; T _{cp} -2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)