

Product datasheet for **GA203213**

Mouse *Pik3c2a* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Pik3c2a</i> CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-4-phosphate 3-kinase catalytic subunit type 2 alpha
Symbol:	<i>Pik3c2a</i>
UniProt ID:	Q61194
Synonyms:	Cpk-m; PI3KC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)